

Communication Protocol

Introduction

This document provides the information needed to communicate with the OEM Valve Hub board through direct UART communication.

Document status & Revision history

Version	Author	Release date	Comments
v1.0.0	Camille Malek	05/04/2024	launch version

Serial connection settings

Baud rate: 230400

Data bits: 8

Stop bit: 1

Parity: none

Termination character: '\n'

Syntax

Command syntax

char 0: '<' to start the query

char 1 to 5: command name

char 6: '?' to read, '!' to write

then ':' to start a value. Can iterate over many arguments

Error handling

In an answer, after the read/write character, '[xx]' with xx 2 hexadecimal numbers are also sent and indicating the error code associated with the request. '00' means non error. The following error codes are:

Error code	Meaning
00	No error
C0	Channel error: wrong channel requested
L0	Locking error: you do not have writing access to this parameter
I0	Impossible command: this query can not be processed
P0	Pause error: this command can not be processed while pause is set to 1
U0	Command incompatible with universal sensor connected to this channel
NU	Command incompatible with non universal (=classic) sensor connected to this channel
B0	Argument value out of bound

List of commands

Parameter	Mandatory arguments	Arguments	W	R	Number of characters returned	Example query	Typical answer	Note
VALVE	int : channel (1 to 4)	bool : valve state	X	X	17	<VALVE?:4 <VALVE!:4:1	>VALVE? 00 04:01 >VALVE! 00 04:01	
IDN		str : device name		X	22	<_IDN_?	>_IDN_? 00 OEMVALVES_	
DEVSN		str : SN		X	18	<DEVSN?	>DEVSN? 00 48V111	
FIRMV		str : firmware version		X	21	<FIRMV?	>FIRMV? 00 v01.03.01	
RESET						<RESET		reset firmware
VALVS		int : 16 bits register for all valve states(0 to 65535)	X	X	17	<VALVS? <VALVS!:65535	>VALVS? 00 65535 >VALVS! 00 65535	
PINGA				X	17	<PINGA?	>PINGA? 00 65535	
PAUSE		int : pause status	X	X	14	<PAUSE?	>PAUSE? 00 00	Pause in ESI

						<PAUSE!:1	>PAUSE! 00 01	0 by default keep state of activated valves
STOP_		int: stop status	X	X	14	<STOP_? <STOP_!:1	>STOP_? 00 00 >STOP_! 00 01	Stop in ESI 0 by default force all valve state to 0